

Development and Implementation of a Novel 3D Culture Technology

Sarah-Louise Ryan^{1,2}, Anthony M. Davies^{1,3}

¹Vale Life Sciences Pty., Translational Research Institute, Brisbane, Australia.

²Cancer and Ageing Research Program, Queensland University of Technology, Brisbane, Australia.

³Irish National Centre for High Content Screening and Analysis, Institute of Molecular Medicine, Trinity College Dublin, Ireland.

Introduction

Specific limitations associated with 2D cell culture have been discerned; such as the loss of tissue-specific architecture, mechanical and biochemical cues, and cell to cell interactions^{1,2}. To address these issues, many are now turning to 3D-based cellular assay systems that permit the formation of multicellular structures (MCS). It has been shown that in contrast to traditional 2D methods, cells cultured in a 3D format may exhibit unique bio-chemical and morphological features similar to their corresponding tissues *in vivo*³. Disease models possessing increased clinical relevance may inform and improve therapeutic intervention. We at Vale Life Sciences, based at the Translational Research Institute, have developed a 3D culture technology that is cost effective, scalable, rapidly deployable, user-friendly and compatible with imaging and liquid handling systems, as well as flow cytometry^{4,5}. Happy Cell Advanced Suspension Medium® (HC) is a low viscosity liquid matrix which promotes the formation of 3D MCS by permitting natural cellular migration, aggregation and proliferation, without restrictions imposed by artificial matrices⁴. HC has also been shown to facilitate the expansion and growth of single cells and enhance the viability of primary cells⁶. When no longer required this suspension system can be inactivated, thereby simplifying cell harvesting and downstream processing. Here, we demonstrate the ability of HC to generate *in vitro* 3D MCS which are compatible with a number of downstream processes.

Methods



HC is supplied as a 4X concentrate and is readily diluted to the required working concentration with standard cell culture medium and desired additives. A 1X concentration is appropriate for most applications. HC is available in a range of culture medium bases and has been designed to sustain cells in long term culture. To use, simply collect cells by centrifugation, re-suspend in HC prepared to the required working concentration, and plate out in low adhesion culture vessels. Using a variety of cell types we have demonstrated its suspension properties, its promotion of cellular viability, and its ability to facilitate the generation of both cellular aggregates and clonal spheroids. Designed to work in conjunction with HC, the complementary Inactivation Solution (IS) facilitates sedimentation and recovery of 3D MCS for downstream processes such as imaging or protein and nucleic acid investigation.

Data




Figure 1: 3D structure formation in HC over a period of 5 days. (IN Cell Analyzer, 4X magnification).




Figure 2: (A) LNCaP and (B) RWPE-1 cells cultured in HC. (Cytell Imaging System, 4X magnification).





Figure 3: 7 day viability and proliferation assessment of Jurkat cells cultured simultaneously in HC and complete culture medium. Cellular viability was significantly increased in cells cultured in HC versus those cultured in standard culture medium (A). A significant increase in HC culture proliferation was also observed after a 7 day period ($p < 0.05$) (B). Data graphed as mean \pm SEM (n=3). Statistical analysis was performed using an unpaired two tailed Welch's t test ($p < 0.05$, $**p < 0.01$).





Figure 4: H460 cell count at surface of well after addition of IS. Mean \pm SEM (n=3). A: Unpaired two tailed Welch's t test ($\delta p < 0.0001$ HC plus IS versus complete medium; $\epsilon p < 0.001$ HC plus IS versus complete medium; $\phi p < 0.01$ HC plus IS versus complete medium; $\psi p < 0.05$ HC plus IS versus complete medium). B: Paired two tailed Welch's t test ($*p < 0.05$; $**p < 0.01$ vs. 0 min).




Figure 5: Presence of necrotic core in larger H460 3D MCS. Cytquant viable cell stain, approx. 150 μ m in diameter (A & B). Hoechst nuclear stain and PI dead cell identifier, approx. 300 μ m in diameter (C & D). (Cytell Imaging System, 10X magnification).




Figure 6: Number of proteins in common and unique to each sample.

Figure 7: Heat map representation of protein expression in 2D and 3D samples. 1337 proteins represented. Scale: 10⁶ based on total Intensity Chi² analysis.

Table 1: 2D upregulated pathways	Table 2: 3D upregulated pathways
Translation	Oxidation reduction
Cell cycle	Apoptosis
RNA processing	Macromolecular complex assembly
Protein localisation	Proteolysis
Protein transport	Generation of precursor metabolites and energy
Proteolysis	Macromolecular catabolic process
Macromolecule catabolic process	Autophagy
Macromolecule complex assembly	Intracellular transport
Regulation of cellular protein metabolic process	Response to organic substance
mRNA processing	Modification-dependant protein catabolic process
RNA splicing	
Cellular response to stress	
Cellular macromolecule localisation	

Conclusion

Happy Cell ASM® has been tested and utilised by major pharmaceutical and biotech companies, contract research organisations and academic institutions. It is featured in an ever increasing number of research communications and is a viable alternative to solid gel matrices, scaffolds, micro-patterned surfaces and hanging drop systems for 3D cell culture. For more information please visit www.valelifesciences.com.

Vale Life Sciences Pty., Level 3, Translational Research Institute, Woollongabba, QLD 4102.
Corresponding email address: enquiries@valelifesciences.com

References

- Lee, J. M. et al. (2013). A three-dimensional microenvironment alters protein expression and chemosensitivity of epithelial ovarian cancer cells *in vitro*. *Lab Invest*, 93(5), 528-542.
- Nelson, C. M., & Bissell, M. J. (2006). Of extracellular matrix, scaffolds, and signaling: tissue architecture regulates development, homeostasis, and cancer. *Annu Rev Cell Dev Biol*, 22, 287-309.
- Kunz-Schughart, L. A. et al. (1998). Multicellular Spheroids: a three-dimensional *in vitro* culture system to study tumour biology. *International Journal of Experimental Pathology*, 1-23.
- Ryan, S. L. et al. (2016). Drug Discovery Approaches Utilizing Three-Dimensional Cell Culture. *Assay Drug Dev Technol*, 14(1), 19-28.
- Takahashi, Y. et al. (2015). A Refined Culture System for Human Induced Pluripotent Stem Cell-Derived Intestinal Epithelial Organoids. *Stem Cell Reports*, 10(1), 314-328.
- Kulasingham, A. et al. (2016). Short term ex-vivo expansion of circulating head and neck tumour cells. *Oncotarget*, 7(37), 60101-60109.