



Happy Cell[®] ASM Inactivation Solution

Product Information Sheet

About Happy Cell[®] ASM Inactivation Solution

Happy Cell[®] ASM Inactivation Solution is a liquid reagent used in conjunction with Happy Cell[®] ASM. It facilitates the recovery of 3D multi-cellular structures by irreversibly disrupting the Happy Cell[®] ASM suspension polymer complex.

Catalogue Number: VHCIS

Components

- 1 x 2 mL vial of Inactivation Solution (5 mg/mL)

Additional Items Required

- Micropipettes
- A 37°C cell culture incubator

Storage and Expiry

Stable until expiry date on bottle if stored at 2-8°C. DO NOT FREEZE.

Ship at ambient temperature.

Safety warnings and precautions

For research use only.

Not recommended or intended for diagnosis of disease in humans or animals.

Do not use internally or externally in humans or animals.

All chemicals should be considered as potentially hazardous. We therefore recommend this product be handled only by persons trained in laboratory techniques and that it is used in accordance with



the principles of good laboratory practice. Wear suitable protective clothing such as laboratory overalls, safety glasses and gloves. Care should be taken to avoid contact with skin or eyes. In the case of contact with skin or eyes wash immediately with water. See material safety data sheet(s) and/or safety statement(s) for specific advice.